

PLACENTAL TISSUE MATRIX SURGICAL TECHNIQUE GUIDE



Storage & Handling



Product can be stored at ambient temperature until the expiry date.

+



CRYO PRODUCT is shipped in a cryo-cooler. Shipping box must be kept room temperature & sealed until use.

+



CRYO PRODUCT container has a 72 hour shelf life. Contact customer service if longer storage time is required.

+



CRYO PRODUCT can be stored at -20°C for up to 6 months, or -40°C until the expiry date.

Preparation & Application

+



CRYO PRODUCT ONLY remove product package from cold storage 15 min prior to use to thaw.

1.



Using aseptic technique, peel open clear product pouch & present vial onto sterile field.

2.



If mixing is desired*, remove cap & add sterile saline or other fluid, and replace cap.

3.



Shake vial to mix completely.

4.



Open vial & draw into syringe using 18G needle (or similar). Change to 21-23G to inject.

5.



Inject at a 45° angle, clockwise throughout the site to provide maximum coverage.

6.



Inject slowly to allow graft to settle into the patient's connective tissues.

7.



Apply the graft directly into the area of damaged or inadequate connective tissue.

* Mix at up to a 2:1 ratio, depending on the preferred total injectable volume for the intended site. Product may be combined with sterile saline or other fluid. Larger volumes may need to be mixed in the syringe.

Dosing


Post-Treatment



Dosing is determined by the physician, depending on the amount or extent of damaged or inadequate connective tissue.



After use, affix Tissue ID Label from box to the included Tracking Card, and return to Skye customer service.



Follow standard post-treatment care protocol.



Physicians may prescribe anti-inflammatories or other medications as they see fit.